

Enhance your hospitals Perioperative efficiency with this Output Adapter/Multiplier

Use with select monitors to simultaneously provide one BP and three ECG analog high level signals on phone Jacks.

Input from
the select MONITOR

order


ME590257F for GE 7 pin mini DIN
ME590257P for GE PDM
ME590257L for GE/Marquette 6 pin LEMO
ME590257 for Datex AS/3 M-PT 4 pin LEMO
ME590257H for Siemens SC7000
SC9000XL or Draeger Delta/Delta XL

ULTRASOUND ECG

ME590132F for GE Vivid i
ME590032R for GE Vivid 5 or 7
ME590151R for GE Vivid 3
ME590338J for Philips/HP Sonos
ME590335 for Siemens
Call us for a Toshiba cable

DEFIBRILLATOR ECG

ME590006K for LP12, LP15 or LP20
ME590335 for 6 pin AAMI (Zoll)
ME590010X for Philips/HP 12 pin



BALLOON PUMP ECG & BP

These two cables (12" long) go to a small junction box that combines them into one 14 foot cable to run High level ECG and IBP from this output to the auxiliary inputs of a Balloon Pump. Both ends of the cable are identical.

Cable Part Number ME590265P



Designed, Manufactured and Distributed by
MAGUIRE Enterprises, Inc.
800-548-9686
www.MaguireEnterprises.com

The ME590257 series does not include the phone plug cables shown here. They must be ordered separately and specifically for your application.